PACKAGING

Ref.: 101-0467	Cont.: 4 x 1 mL
----------------	-----------------

Store at 2 - 8° C

PRODUCT CHARACTERISTICS

The control is a lyophilised human serum used to evaluate the precision and accuracy of ASO, CRP and RF determinations in human serum by turbidimetry method.

REAGENTS

ASO/CRP/RF Control H	Human serum. With a high concentration of ASO, CRP and RF.
----------------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Control have been standardized against the Reference Material ERM-DA 472/IFCC for CRP, International Reference Standard from NIBSC 64/002 for RF and International Standard from NIBSC 97/662 for ASO.

PREPARATION

Reconstitute accurately with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 10 days at 2-8° C or 6 months at -20° C.

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Store at 2 - 8° C

PRODUCT CHARACTERISTICS

The control is a lyophilised human serum used to evaluate the precision and accuracy of ASO, CRP and RF determinations in human serum by turbidimetry method.

REAGENTS

ASO/CRP/RF Control H	Human serum. With a high concentration of ASO, CRP and RF.
----------------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Control have been standardized against the Reference Material ERM-DA 472/IFCC for CRP, Rheumatoid Arthritis Serum WHO 64/1 for RF and ASO 1st International Calibrator WHO for ASO.

PREPARATION

Reconstitute accurately with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 10 days at 2-8° C or 6 months at -20° C.

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