**CRP/PCR-Ultra Control**  
Turbidimetry / Turbidimetria

*Quantitative determination of CRP-ultra*

**PACKAGING**

Ref.: 101-0510	Cont.: 1 x 2 mL
----------------	-----------------

Store at 2 - 8° C

**PRODUCT CHARACTERISTICS**

The control is a human serum containing human CRP used to evaluate the precision and accuracy of CRP-ultra determinations in human serum by turbidimetry method.

**REAGENTS**

U-CRP Control	Human serum. With a low concentration of CRP. Preservative. The CRP concentration is stated on the vial label.
---------------	--

**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**CALIBRATION**

The concentration value of the control has been standardized against the Reference Material ERM-DA472/IFCC

**PREPARATION**

Ready for use.

**STORAGE AND STABILITY**

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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**CRP/PCR-Ultra Control**  
Turbidimetry / Turbidimetria

*Quantitative determination of CRP-ultra*

**PACKAGING**

Ref.: 101-0510	Cont.: 1 x 2 mL
----------------	-----------------

Store at 2 - 8° C

**PRODUCT CHARACTERISTICS**

The control is a human serum containing human CRP used to evaluate the precision and accuracy of CRP-ultra determinations in human serum by turbidimetry method

**REAGENTS**

U-CRP Control	Human serum. With a low concentration of CRP. Preservative. The CRP concentration is stated on the vial label.
---------------	--

**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**CALIBRATION**

The concentration value of the control has been standardized against the Reference Material ERM-DA472/IFCC

**PREPARATION**

Ready for use.

**STORAGE AND STABILITY**

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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