IgE Control
Turbidimetry / Turbidimetria

IgE Control
Turbidimetry / Turbidimetria

Quantitative determination of IgE

PACKAGING

Ref.: 101-0480	Cont.: 1 x 2 mL
----------------	-----------------

Store at 2 - 8° C.

PRODUCT CHARACTERISTICS

The control is a liquid IgE solution from human origin containing IgE used to evaluate the precision and accuracy of IgE determinations in human serum by turbidimetry method.

REAGENTS

IgE- Control	Human IgE serum. With a low concentration level of IgE. The IgE concentration is stated on the vial label.
--------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against the International Reference Preparation of Human Serum IgE 75/502 (2nd IRP, 1981).

PREPARATION

Ready to use

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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Quantitative determination of IgE

PACKAGING

Ref.: 101-0480	Cont.: 1 x 2 mL
----------------	-----------------

Store at 2 - 8° C.

PRODUCT CHARACTERISTICS

The control is a liquid IgE solution from human origin containing IgE used to evaluate the precision and accuracy of IgE determinations in human serum by turbidimetry method.

REAGENTS

IgE -Control	Human IgE serum. With a low concentration level of IgE. The IgE concentration is stated on the vial label.
--------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against the International Reference Preparation of Human Serum IgE 75/502 (2nd IRP, 1981).

PREPARATION

Ready to use

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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