**Microalbumin Control /  
Microalbúmina Control**  
Turbidimetry / Turbidimetria

*Quantitative determination of microalbumin*

**PACKAGING**

|                |                 |
|----------------|-----------------|
| Ref.: 101-0483 | Cont.: 1 x 2 mL |
|----------------|-----------------|

Store at 2 - 8° C

**PRODUCT CHARACTERISTICS**

The control is a solution containing human albumin used to evaluate the precision and accuracy of microalbumin determinations in human urine by turbidimetry method.

**REAGENTS**

|                  |   |
|------------------|---|
| μALB-<br>Control | Human albumin solution. With a low concentration level of albumin. Preservative. The albumin concentration is stated on the vial label. |
|------------------|---|

**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**CALIBRATION**

The concentration value of the Control has been standardized against the International Reference Material CRM 470/RPPHS (BCR).

**PREPARATION**

Ready for use

**STORAGE AND STABILITY**

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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Store at 2 - 8° C

**PRODUCT CHARACTERISTICS**

The control is a solution containing human albumin used to evaluate the precision and accuracy of microalbumin determinations in human urine by turbidimetry method.

**REAGENTS**

|                  |  |
|------------------|--|
| μALB-<br>Control | Human albumin solution. With a low concentration level of albumin. Preservative. The albumin concentration is stated on the vial label |
|------------------|--|

**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**CALIBRATION**

The concentration value of the Control has been standardized against the International Reference Material CRM 470/RPPHS (BCR).

**PREPARATION**

Ready for use

**STORAGE AND STABILITY**

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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