

PACKAGING

Ref.: 101-0366	Cont.: 12 x 50 mL
Ref.: 101-0391	Cont.: 4 x 50 mL

Store at 2 - 8° C

CLINICAL SIGNIFICANCE


Glucose is a major source of energy for most cells of the body; insulin facilitates glucose entry into the cells.

Diabetes is a disease manifested by hyperglycemia; patients with diabetes demonstrate an inability to produce insulin^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

PRINCIPLE OF THE METHOD

Hexokinase (HK) catalyzes the phosphorylation of glucose to glucose-6-phosphate (G6P) by ATP. The formed glucose-6-phosphate is reduced to 6-phosphogluconate in the presence of glucose-6-phosphate dehydrogenase (G6P-DH) with the subsequent reduction of NAD to NADH:



The increase in concentration of NADH is proportional to the glucose concentration in the sample^{1,2}.

REAGENTS

R 1 Buffer	TRIS pH 7.5	4 mmol/L
	ATP	2.1 mmol/L
	Mg ²⁺	0.8 mmol/L
R 2 Enzymes	NAD ⁺	2 mmol/L
	Hexokinase (HK)	1000 U/L
	Glucose-6-phosphate (G6P-DH)	1000 U/L
GLUCOSE CAL	Glucose aqueous primary standard 100 mg/dL	

Optional (not included in the kit)

Contro-N	Ref.: 101-0252	4 x 5 mL	Lyophilized human control serum
	Ref.: 101-0083	20 x 5 mL	
Contro-P	Ref.: 101-0253	4 x 5 mL	Lyophilized human control serum
	Ref.: 101-0084	20 x 5 mL	

PREPARATION

Working reagent (WR): Dissolve (→) the contents of one vial of R 2 Enzymes in one bottle of R 1 Buffer.

Cap and mix gently to dissolve contents.

The reagent is stable after reconstitution 1 month in the refrigerator (2 - 8° C) or 7 days at room temperature (15 - 25° C).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2 - 8° C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm \geq 0.30.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

- Serum or plasma free of hemolysis¹.

- Serum should be removed from the clot as quickly as possible.

- Urine¹.

Stability of the sample: Glucose in serum or plasma is stable at 2 - 8° C for 3 days.

PROCEDURE

Notes: CHRONOLAB SYSTEMS has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

GLUCOSE CAL: Proceed carefully with this product because due its nature it can get contaminated easily.

Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.

Use clean disposable pipette tips for its dispensation.

1. Assay conditions:
Wavelength: 340 nm
Cuvette: 1 cm. light path
Temperature: 37° C / 15 - 25° C
2. Adjust the instrument to zero with distilled water.
3. Pipette into a cuvette:

	Blank	Standard	Sample
WR (mL)	1.0	1.0	1.0
Standard ^(Note 1-2) (µL)	--	10	--
Sample (µL)	--	--	10

4. Mix and incubate for 5 min. at 37° C or 10 min. at room temperature (15 - 25° C).
5. Read the absorbance (A) of the samples and calibrator, against the Blank.

CALCULATIONS

$$\frac{(A)Sample - (A)Blank}{(A)Standard - (A)Blank} \times 100 \text{ (Calibrator conc.)} = \text{mg/dL glucose in the sample}$$

Conversion factor: mg/dL x 0.0555 = mmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures.

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

Serum or plasma:

$$60 - 110 \text{ mg/dL} \quad \cong \quad 3.33 - 6.10 \text{ mmol/L}$$

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0.16 mg/dL to *linearity limit* of 600 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl (9 g/L) and multiply the result by 2.

Precision:

Mean (mg/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
	99.5	244	98.0	247
SD	0.83	1.70	1.60	3.75
CV (%)	0.83	0.70	1.63	1.51

Sensitivity: 1 mg/dL = 0.0036 A.

Accuracy: Results obtained using CHRONOLAB reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0.9891

Regression equation: $y=9877x + 0.9366$

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Hemoglobin up to 19 g/L and bilirubin up to 100 mg/L, do not interfere¹.

A list of drugs and other interfering substances with glucose determination has been reported by Young et. al^{3,4}.

BIBLIOGRAPHY

1. Kaplan L.A. Glucose. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1032-1036.
2. Trinder P. Ann Clin Biochem 1969; 6: 24-33.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.