**Lp(a) Calibrator /
Lp(a) Calibrador**
Turbidimetry / Turbidimetria

Quantitative determination of Lipoprotein (a) (Lp(a))

PACKAGING

Ref.: 101-0531	Cont.: 1 x 1 mL
----------------	-----------------

Store at 2 - 8° C.

PRODUCT CHARACTERISTICS

The calibrator is a lyophilised serum from human origin containing a Lp(a) concentration suitable to calibrate Lp(a) determinations by turbidimetry method.

REAGENTS

Lp(a)-CAL	Calibrator. Human serum. Sodium azide 0.95 g/L. Lp(a) concentration is stated on the vial label.
-----------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Calibrator has been standardized against an Internal Reference Material.
Reference Method: ELISA method stated in *J. Lipid Research 30, 651-662, 1989*.

PREPARATION

Reconstitute with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring to room temperature for about 10 minutes before use.

STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 1 month at 2-8° C or 3 months at -20° C.

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