HbA_{1c} Multicalibrador / Multicalibrador

Latex turbidimetry / Turbidimetria Látex
Quantitative determination of glycated hemoglobin (HbA_{1c})

PACKAGING

Ref.: 101-0559	Cont.: 6 x 2 mL
----------------	-----------------

Store 2 - 8° C.

PRODUCT CHARACTERISTICS

The HbA_{1c} CAL is a multicalibrator from human origin containing different HbA_{1c} concentrations suitable to calibrate HbA_{1c} determinations by turbidimetry method.

REAGENTS

HbA _{1c} Calibrator Level 1	Alkaline Solution, NaOH 0.4 %
HbA _{1c} Calibrator Level 2	Human HbA _{1c} . Contains preservatives. The concentration in g/dL is indicated on the insert.
HbA _{1c} Calibrator Level 3	
HbA _{1c} Calibrator Level 4	
HbA _{1c} Calibrator Level 5	
HbA _{1c} Calibrator Level 6	

PRECAUTIONS

Calibrator Level 1: Irritant. Contains sodium hydroxide, and may cause severe burns. Use gloves. In the event of accident, flush the effected area with plenty of water and seek immediate medical attention.

Calibrator Level 2 to 6: Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

Calibrator values have been standardized against an Internal Reference Material HbA_{1c} traceable to the HPLC method for HbA_{1c} and Drabkin's method for Total Hb.

PREPARATION

Reagents: Ready to use. Mix well by inversion to ensure homogeneity.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2 - 8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

PROCEDURE

Calibrators do not require pretreatment.

Total Hb test: This test is performed with HbA_{1c} Calibrator Level 1.

HbA_{1c} test: This test is performed with HbA_{1c} Calibrator Level 1 to 6.

* This value sheet is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C, etc.) following the lot number.

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