### PACKAGING

Ref.: 101-0560	Cont.: 2 x 0.5 mL
----------------	-------------------

Store 2 - 8° C.

### PRODUCT CHARACTERISTICS

The HbA<sub>1c</sub> CONTROL 1/2 is multicontrol from human origin containing two different HbA<sub>1c</sub> concentrations suitable to control the precision and accuracy of HbA<sub>1c</sub> determinations by turbidimetry method.

### REAGENTS

<b>HbA<sub>1c</sub> Control Level 1</b>	Human blood. Contains preservatives. The HbA <sub>1c</sub> concentration in g/dL is indicated on the insert.
<b>HbA<sub>1c</sub> Control Level 2</b>	

### PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

### CALIBRATION

Control values have been standardized against an Internal Reference Material HbA<sub>1c</sub> traceable to the HPLC method for HbA<sub>1c</sub> and Drabkin's method for Total Hb.

### PREPARATION

1. Remove the cap from the control vials.
2. Reconstitute → with 0.5 mL of double deionised water.
3. Cap the vials and mix thoroughly, avoiding foaming.
4. Bring at room temperature for at least 15 minutes.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2 - 8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Stability of the reconstituted controls:** 1 month when stored at 2 - 8° C. *Do not freeze the reconstituted controls.*

### PROCEDURE

**Controls are treated the same as samples.**

Mix into separate tubes, 10 µL of each reconstituted control with 400 µL of Haemolizing Reagent (R3) ( 1/41 Dilution).

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