## Serum control multicomponent for clinical chemistry assays IVD

### PACKAGING

Contro-N	
Ref.: 101-0252	Cont.: 4 x 5 mL
Ref.: 101-0083	Cont.: 20 x 5 mL
Contro-P	
Ref.: 101-0253	Cont.: 4 x 5 mL
Ref.: 101-0084	Cont.: 20 x 5 mL

Store at 2 - 8° C.

### PRODUCT CHARACTERISTICS

Contro-N or Contro-P is a human lyophilised serum. With most constituent concentrations and activities in the normal or pathologic range. It's intended for control of accuracy for use with manual and automated analytical procedures.

### REAGENTS

Human serum. Enzymes (human, bovine and porcine). Nonprotein constituents and bacteriostatic agents. The concentration / activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

### PRECAUTIONS

The individual units comprising lot or lots of components from human origin, have been tested for HBsAg and antibodies to HIV and HCV found negative / non reactive by an FDA approved test. However, handle cautiously as potentially infectious.

### PREPARATION

-Using a Class A volumetric pipette, carefully add the exact amount of distilled or deionized water as specified on the vial label to the control material.  
-Swirl the vial gently and allow to stand at room temperature for 5 to 10 minutes.  
-Swirl contents until homogeneous. DO NOT VORTEX. Remove the sample required for testing. Improper handling and/or storage of the control can affect results. Inaccurate reconstitution of the control and errors in assay technique can cause erroneous results.

### STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.  
Store tightly capped when not use.  
- After reconstitution is stable for:  
At 15° C to 25° C: 12 hours (\* 8 hours Bilirubin Total; 4 hours Bilirubin Direct, ACP and ALP)  
At 2-8° C: 5 days (\*24 hours Bilirubin Total, Bilirubin Direct, ACP and ALT).  
At -15° C to -25° C: 1 month (\*2 weeks Bilirubin Total, Bilirubin Direct, ACP and ALP)

### BIBLIOGRAPHY

1. Council Directive (2000/54EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.
2. International Federation of Clinical Chemistry (IFCC) Education Division, Expert Panel of Quantities and Units: *A Protocol for the Conversion of Clinical Laboratory data*, Journal of Automatic Chemistry Vol. 11, No 5 (Sept – Oct 1989), pp. 223-226

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