PACKAGING

Ref.: 101-0774	Cont.: 4 x 0.5 mL (2 levels)
----------------	------------------------------

Store 2 - 8° C.

INTENDED USE

This product is for the purpose of monitoring accuracy and precision in the quantitative determination of human hemoglobin A1c (HbA_{1c}) in blood by automated immunoassay.

PRODUCT CHARACTERISTICS

The lyophilized HbA_{1c} controls are hemolysates prepared from packed human erythrocytes. The controls provide two levels of HbA_{1c}, one level in the normal range and the other level in the elevated range. Stabilizers are added to maintain hemoglobin in the reduced state providing complete control of the HbA_{1c} procedure.

REAGENTS

HbA_{1c} Control Level 1	Human blood. Contains preservatives. The HbA _{1c} concentration is indicated on the vial label.
HbA_{1c} Control Level 2	

PREPARATION

Reconstitute vials with 0.5mL deionized water. Gently mix for 10 minutes. Observe for undissolved material. The reconstituted controls may be dispensed in 0.1 mL aliquots, sealed tightly and frozen at -20°C.

The lyophilized HbA_{1c} controls should be assayed in the same manner as blood specimens including the hemolysate procedure.

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious. Avoid contact with skin and mucous membranes.

MATERIAL REQUIRED AND NOT PROVIDED

Hemoglobin A1c Reagent kit, from Chronolab Ref: 101-0771
Hemoglobin A1c Calibrator kit, from Chronolab Ref: 101-0773

STORAGE AND STABILITY

1. Store controls at 2-8° C. Stable until expiration date if sealed tightly. **PROTECT FROM LIGHT AND HEAT.**
2. Reconstituted controls retain their assigned values for at least three months if frozen. If not frozen, the reconstituted controls are stable at least one month if stored at 2-8° C and sealed tightly.
3. Do not freeze and thaw more than once.

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