**Microalbumin Control /
Microalbúmina Control**
Turbidimetry / Turbidimetria

**Microalbumin Control /
Microalbúmina Control**
Turbidimetry / Turbidimetria

Quantitative determination of microalbumin

PACKAGING

Ref.: 101-0483	Cont.: 1 x 2 mL
----------------	-----------------

Store at 2 - 8° C

PRODUCT CHARACTERISTICS

The control is a solution containing human albumin used to evaluate the precision and accuracy of microalbumin determinations in human urine by turbidimetry method.

REAGENTS

μALB- Control	Human albumin solution. With a low concentration level of albumin. Preservative. The albumin concentration is stated on the vial label.
------------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against the International Reference Material ERM-DA 470K/IFCC (BCR).

PREPARATION

Ready for use

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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Store at 2 - 8° C

PRODUCT CHARACTERISTICS

The control is a solution containing human albumin used to evaluate the precision and accuracy of microalbumin determinations in human urine by turbidimetry method.

REAGENTS

μALB- Control	Human albumin solution. With a low concentration level of albumin. Preservative. The albumin concentration is stated on the vial label
------------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Control has been standardized against the International Reference Material ERM-DA 470K/IFCC (BCR).

PREPARATION

Ready for use

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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