Quantitative determination of Apolipoproteins (APO)

PACKAGING

Ref.: 101-0503	Cont.: 1 x 1 mL
----------------	-----------------

Store at 2 - 8° C.

PRODUCT CHARACTERISTICS

APO CONTROL is a lyophilised serum from human origin containing apolipoprotein concentrations suitable to evaluate the accuracy and precision of apolipoproteins determinations by turbidimetry method.

REAGENTS

APO-CONTROL	Control. Human serum. Sodium azide 0.95 g/L. APO concentrations are stated on the list below.
-------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Control are traceable to Secondary Reference Material Apo A-I WHO / IFCC (SP1-01) and Apo B WHO / IFCC (SP3-07), and Internal References Materials for APO C-II, C-III and E.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming.

The resulting control solution should be used immediately or stored at 2-8° C.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8° C and contaminations are prevented during their use. Do not use reagents over the expiration date.

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This Control data sheet for lot is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number. This product was manufactured, packaged and tested in accordance to CHRONOLAB, S.L. Quality System and meets all product specifications. Our Quality System has been approved by an accredited registering body to ISO 9001.

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